

M SERIES HUMIDIFICATION PRODUCT REFERENCE GUIDE

This reference guide provides identification and ordering information for M Series humidifiers and water chambers.

HUMIDIFIER AND CHAMBER COMPATIBILITY



All humidifiers with an underlined 5 are compatible with J-tube chamber part number 1003758.

Humidifiers available to order:

- 1049109 DC Standalone Heated Humidifier
- 1051081 DC Standalone Passover Humidifier
- All integrated device packages

Humidifiers that are obsolete but still compatible:

- 1037652 AC Standalone Heated Humidifier



J-tube chamber part number 1003758.



All humidifiers with the 5 not underlined are compatible with side-by-side chamber part number 1003757.

Humidifiers that are obsolete but still compatible:

- 1037652 AC Standalone Heated Humidifier
- 1029090 AC Standalone Passover Humidifier



Side-by-side chamber part number 1003757.

CHAMBER KIT



The universal kit provides one chamber bottom with both a J-tube and side-by-side chamber top.

- 1003756 Universal RP Tank Kit

ADDITIONAL ACCESSORIES



Inlet/outlet seal for J-tube.

Replacement seal for J-tube chamber

- 1003784 Inlet/outlet seal for J-tube, 2 per pack

Replacement seal for side-by-side chamber

- 1028116 Inlet/outlet seal for side-by-side, 1 per pack
- 1029357 Inlet/outlet seal for side-by-side, 10 per pack

Chamber O-ring, part number 1029358 is obsolete and should be replaced with new chamber part number 1003757.



Inlet/outlet seal for side-by-side.