



Children's  
Medical  
Ventures

A NOVAMETRIX COMPANY

# Specialty Feeding Products

Imported from Japan



## Cleft Palate Nipple

- Specially designed nipple for feeding of cleft palate and surgical patients, babies with other oral abnormalities, or infants with poor sucking strength.
- Cross-cut nipple has a thick and a thin side for easier sucking
- One-way valve prevents liquid from flowing back into the bottle and only allows milk to flow when sucked by the baby
- Special air vents help caregiver to adjust the rate of flow
- May be purchased separately or with baby bottle.
- The nipple and one-way valve fit other baby bottles, but may have some leaking or nipple collapse
- Latex free!

## Baby Bottle

- Specifically made and shaped to fit the air vents and one-way valve of the cleft palate nipple

## Food Feeder

- 120ml squeezable bottle with an attached spoon
- Assists in the weaning process from bottle to spoon or to aid problem feeders.
- Suitable for juice, soups and milk.

275 Longwater Drive  
Norwell, MA 02061  
[www.childmed.com](http://www.childmed.com)

# Specialty Feeding Products

Imported from Japan

## Specifications

### Cleft Palate Nipple:

- Fits standard baby bottle or volu-feeder, but works best with Baby Bottle
- Material: non-latex isoprene nipple, polypropylene valve
- Individually packaged
- Semi-transparent "caramel" color
- Sanitize by placing in boiling water for two minutes
- Do not autoclave or steam sterilize

### Baby Bottle:

- Material (bottle, ring & cap): Polypropylene, non-latex
- Sanitize by placing in boiling water for five minutes

### Food Feeder Bottle:

- Material: Polypropylene, non-latex
- Sanitize by placing in boiling water for five minutes.

## Instructions for Use

- A booklet is available that explains the sucking movement of a baby and the design of the Pigeon Specialty Feeding products.
- Please read "Instructions for Use" prior to using.

Children's Medical Ventures is an authorized distributor of Pigeon Specialty Feeding Products. Pigeon Corporation, based in Tokyo, Japan, is a leading manufacturer of baby care products and is committed to high quality innovative products and research.

## Ordering Information

ITEM NO.	DESCRIPTION	QTY.
98072	Food Feeder Bottle	12/CS
98073	Cleft Palate Nipple	20/CS
98074	Baby Bottle w/ two (2) cleft palate nipples	10/CS

**TO ORDER**  
**CALL 800-377-3449**  
**or 781-871-6226**  
**or FAX 203-294-2330**

[www.childmed.com](http://www.childmed.com)

**275 Longwater Drive**  
**Norwell, MA 02061**