

# Cozycare™ Bunting

## Infant Thermoregulation Wrap



Children's  
Medical  
Ventures

RESPIRONICS®

Children's Medical Ventures' Cozycare™ Buntings are used for providing thermoregulation support and containment in the delivery room, nursery, during transport, during feeding or at home. The unique design combines a layer of DuPont™ Thermolite® material with refastenable closures, making the bunting versatile for many uses.

Cozycare™ Buntings are made from an ultra soft 100% polyester, non-woven material in a unisex print. The material is non-allergenic, non-toxic, flame retardant and the layer of DuPont™ Thermolite® that helps keep the baby warm and comfortable. Cozycare™ Buntings are disposable for hospital use or can be washed several times in a home setting. The buntings come in three sizes and are easily adjustable with soft, quiet refastenable closures.

## Cozycare™ Buntings:

### Provide Thermal Support

- in the delivery room or nursery
- during transport
- during feeding
- at home

### Easy to Use

- disposable in a hospital setting to help reduce linen costs
- washable in a home setting

### Provide Security

- refastenable closures keep baby securely contained even when agitated



### ordering information

item no.	description	quantity
93019-C	Micro Preemie <3 lbs.	25/case
93019-B	Preemie 3–6 lbs.	25/case
93019-A	Newborn 5–10 lbs.	25/case



Customer Service: 1-800-345-6443 or 724-387-4000  
Respironics Europe: +33-1-47-52-30-00  
Respironics Asia Pacific: +81-3-5280-9611

Respironics and the Children's Medical Ventures heart design logo are registered trademarks of Respironics, Inc. and its affiliates. Cozycare is a trademark of Respironics, Inc. and its affiliates. © 2005 Respironics, Inc. and its affiliates. All rights reserved. DuPont and Thermolite are trademarks or registered trademarks of E.I. du Pont de Nemours and Company or its affiliates.