

FileWire for Synergy™

Upload and View Sleep Study Files
from Anywhere in the World



Children's
Medical
Ventures

RESPIRONICS®

Using a secure hosted server, DME's, doctors, clinics and hospitals can now upload and view sleep study files from anywhere in the world. This is made possible through a HIPAA compliant, secure file transfer service from FileWire Technologies.

With a click of a mouse, a sleep lab technician can transmit the results of a patient's sleep study to the File Wire central server. The software finds the various files associated with the study, compresses, encrypts and sends them to the server via the Internet, if available, or a phone line and modem. After the files are uploaded, a message is sent to the recipient indicating that the study is ready for analysis. The study can then be downloaded and scored or interpreted.

FileWire for Synergy™ is:

- Secure - The software provides encryption to ensure data privacy.
- Flexible - Anyone with a computer modem or internet access can use the easy-to-install software.
- Easy - A simple in-box/out-box and send/receive format makes usage a breeze.
- Fast - Even large 30 MB sleep study files can be sent rapidly.
- Economical - No more overnight charges or courier costs. A low monthly fee to upload and download studies.
- Timely - Studies can be processed quickly.

A single click is all that is needed to transfer files.

In Box/Out Box layout makes it easy to understand and operate.

Last Name	First Name	Acq Date	Main Unit	Acq No
DR. VIOLARIS		07/25/2002	1237	28
PENDELTON		11/30/2001	001050	4096
PH STUDY		05/10/2000	001050	10
PH STUDY		05/10/2000	1237	10


**For more information, call FileWire at 301-498-7797
or visit them on the web at www.filewire.net**



Customer Service: 1-800-345-6443 or 724-387-4000
Respironics Europe: +33-1-55-60-19-80
Respironics Asia Pacific: +81-3-5280-9611

Respironics and People. Products. Programs. are registered trademarks of Respironics, Inc. and its affiliates. Synergy is a trademark of Respironics, Inc. and its affiliates.
© 2005 Respironics, Inc. and its affiliates. All rights reserved.

Gyer KB 5/06/05 MCI 4100141 1025857



SLEEP
MANAGEMENT



CHRONIC RESPIRATORY
MANAGEMENT



TOTAL VENTILATION
SOLUTIONS



RESPIRATORY
DRUG DELIVERY



NEONATAL &
INFANT CARE