


Nasal Aspirators

Easy Nasal Suctioning at Your Fingertips

BBG™ Nasal Aspirator (for hospital use)



Bulb-Style Syringe Nasal Aspirator (for hospital and home use)



Children's
Medical
Ventures

RESPIRONICS®

Children's Medical Venture's line of nasal aspirators offer easy and safe suctioning. In the hospital or in the home, the removal of excess mucus or secretions is simple with either the bulb-style syringe or the BBG, which easily connects to standard wall suction. Whichever you choose, you will have **easy nasal suctioning at your fingertips.**

Nasal Aspirators



**BBG™ Nasal Aspirator
(for hospital use)**

Developed in conjunction with Miami Children's Hospital, the BBG™ Nasal Aspirator allows easy and safe suctioning of the nasal or oral cavity. The clear bendable tube allows flexible positioning of the tip for a variety of patients in any position.

Safe Suctioning:

- The bulb end features a graduated, rounded tip to help minimize potential tissue trauma caused by deep suctioning
- Single patient use helps to minimize potential of bacterial colonization

Easy to Use:

- Clear, bendable PVC tubing provides flexibility for suctioning hard to reach places on patients in a variety of positions
- Thumb-control valve allows gentle and controlled suctioning
- Connects to standard wall suction - no special equipment required



**Bulb-Style Syringe Nasal Aspirator
(for hospital and home use)**

This safe, easy to use bulb-style syringe nasal aspirator from The First Years® provides gentle suction to remove excess mucus or oral secretions. The First Years® Aspirator is safer than most bulb syringes: an over-insertion guard helps avoid inserting the tip too far into the baby's nostril and a pull-out plug in the bottom makes it possible to thoroughly clean the interior. This pediatrician-approved nasal aspirator has previously been available commercially to parents, but is now offered to hospitals through Children's Medical Ventures.

Safe Suctioning:

- Features an over-insertion guard to help avoid inserting the tip too far into the nostril
- Single patient use helps to minimize potential of bacterial colonization
- Pediatrician approved

Easy to Use:

- Pull-out plug for thorough, easy cleaning



ordering information

item no.	description	quantity
93032-S	BBG Nasal Aspirator 4.5" Length	50/case
93032-02	BBG Nasal Aspirator 3.25" Length	50/case
01053-1	Bulb Syringe Style Nasal Aspirator	24/case



Customer Service: 1-800-345-6443 or 724-387-4000
 Respiroics Europe: +33-1-55-60-19-80
 Respiroics Asia Pacific: +81-3-5800-0724

BBG is a trademark of Respiroics, Inc. and its affiliates. People. Products. Programs. is a registered trademark of Respiroics, Inc. and its affiliates.
 © 2004 Respiroics, Inc. and its affiliates. All rights reserved.
 The First Years is a trademark of of The First Years, Inc.

1023320 KB 9/21/04

