

Synergy™-E Event Software

For Use with SmartMonitor® 2 & SmartMonitor® 2PS



Children's
Medical
Ventures

RESPIRONICS®

Children's Medical Ventures' Synergy™-E Event Software is a Windows®-based program that simplifies the review of infant apnea data. It provides graphical displays, the ability to define patient information, customizable event labels and a notes section for improved communication.

Our new Synergy-E Windows®-based software is both easy to use and easy to understand. Waveform colors and scaling are adjustable for easy viewing. You can even define patient information, customize event labels and add notes for improved communication. All these features combine to make Synergy-E™ Event Software a powerful and flexible tool for reviewing apnea data.

New Report Design

- Patient demographic data, monitoring parameters and history prints on one page
- 4 concise reports:
 - patient event report
 - equipment report
 - summaries report
 - compliance report

Explanation of Reports

- Patient event report
 - chronological listing of recorded patient events
- Equipment event report
 - displays all monitor settings and parameters
 - chronological listing of all equipment events
- Summaries Report
 - alarm summary
 - event summary
 - heart rate summary
 - ranges of bradycardia events
 - ranges of tachycardia events
 - SpO₂ summary
 - Low SpO₂ % range
 - Low SpO₂ durations
 - associated events - chronological list of events that occur in the same EGRAM
- Compliance report
 - easy-to-read graph of monitor usage, including daily hours of use

Custom Report Headings

- Standard heading data
- User selectable headings
- Customize header/footer options

New Open Text Fields

- Enter virtually unlimited data in text fields:
 - diagnosis
 - medications
 - history
 - comments
 - interpretation
 - recommendations

Synergy-E™ Event Software works with SmartMonitor® 2 and SmartMonitor® 2PS infant apnea monitors.



SmartMonitor® 2PS

ordering information

item no.	description	quantity
1022029	Synergy-E™ Event Software V3.1	1/each
1014557	SmartMonitor® 2PS	1/each
4002	SmartMonitor® 2 w/ PCMCIA Card Slot	1/each
4003	SmartMonitor® 2 w/ PCMCIA Card Slot and Internal Modem	1/each



Customer Service: 1-800-345-6443 or 724-387-4000
 Respironics Europe: +33-1-47-52-30-30
 Respironics Asia Pacific: +81-3-5280-9611

Respironics and the Children's Medical Ventures heart design logo are registered trademarks of Respironics, Inc. and its affiliates. Synergy is a trademark of Respironics, Inc. and its affiliates. © 2005 Respironics, Inc. and its affiliates. All rights reserved.

Gyer KB 10/12/05 MCI 1029764 PN 4100371

