

WeePee® Three Diapers

Designed With The Tiniest Babies In Mind



Medium (for babies 1250-2250 grams)

Small (for babies 800-1250 grams)

Extra Small (for babies <800 grams)



Children's
Medical
Ventures

RESPIRONICS®

Children's Medical Ventures WeePee® Three is the only developmentally appropriate diaper that comes in three sizes. Soft, easy to fasten (and refasten) non-adhesive tabs ensure a comfortable fit. The proprietary dual-function pad is not only super absorbent but also allows urine extraction.

WeePee® Three: the diapers designed with the tiniest babies in mind.

WeePee® Three Diapers

Children's Medical Ventures developed the first diapers specifically for premature infants to help ensure that even the smallest babies enjoy increased comfort and developmentally appropriate positioning. With direction from NICU nurses, WeePee® Three diapers have added even more features:

Three sizes for a custom, comfortable fit for all preemies:

- **Extra Small** (for babies <800 grams)
- **Small** (for babies 800-1250 grams)
- **Medium** (for babies 1250-2250 grams)

Narrow, proprietary dual function pad

- Allows easy urine extraction which eliminates the need for cotton balls or gauze
- Super absorbent to keep moisture away from baby's skin
— Amount absorbed: XS: 25cc; S: 50cc; M: 70cc
- Designed to minimize the risk of hip abduction

Soft, flexible refastenable closures

- Easy to attach and reattach diaper tabs for a custom fit
- Will not stick to baby's skin

Double-stranded, softer elastic leg gathers

- Snugger fit helps minimize leaks
- Softer elastic enhances comfort on delicate skin

Color-coded packaging to help differentiate sizes at a glance

- Purple for extra small
- Teal for small
- Peach for medium



EXTRA SMALL SMALL MEDIUM

Distinct bag colors help to identify the three different WeePee® Three diaper sizes.



Customer Service: 1-800-345-6443 or 724-387-4000
 Respironics Europe: +33-(0)1-55-60-19-80
 Respironics Asia Pacific: +852-234-342-18

ordering information

item no.	description	quantity
1017824	Extra Small (XS) <800 g.	360/case*
1017855	Small (S) 800-1250 g.	360/case*
1017856	Medium (M) 1250-2250 g.	360/case*

*30 diapers per sleeve; 12 sleeves per case

WeePee is a trademark of Respironics Inc and its affiliate.
 ©2004 Respironics, Inc and its affiliates. All rights reserved.

1020146 KB 5/17/04 G

