


# WeeSpecs® Phototherapy Mask

When Baby's Eyes Need Protection


## Preemie

Fits head circumference <28cm




## Small

Fits head circumference 28-34cm



## Medium

Fits head circumference >34cm



Children's Medical Ventures' WeeSpecs® Phototherapy Mask has a unique goggle design that helps to eliminate virtually all light from reaching the baby's eyes during phototherapy or any time light protection is required. **WeeSpecs®...when baby's eyes need protection.**

# WeeSpecs® Phototherapy Mask

WeeSpecs® two-piece design is made from a very soft hypo-allergenic material and can be adjusted for a custom fit. The inside frame of the goggles seals around the skeletal socket, rather than laying flat, reducing pressure on the eyeball which allows a baby's eyes to open and close naturally. WeeSpecs® come in three sizes for a perfect fit.

## Unique Design

- Soft cushioned inner eye “goggles” seal around the skeletal socket rather than laying flat. This reduces pressure on the eyeball and allows baby’s eyes to open and close naturally.



GOGGLE DESIGN

- Outer light blocking material is five layers thick effectively shielding a baby’s eyes during phototherapy or anytime light protection is required.
- Specially designed back-piece conforms to the shape of the baby’s head and soft Velcro® straps keep the mask securely in place.



BACK VIEW

## Three Sizes to Fit Premies and Newborns

- Premie - fits head circumference <28cm
- Small - fits head circumference 28-34cm
- Medium - fits head circumference >34cm

## Specifications

- Flame retardant
- Meets fire hazard classification:
  - FMUSS302 - passed
  - Flame spread - 22
  - California Test 117-passed
- Latex-free



## ordering information

WeeSpecs®

item no.	description	quantity
99071-01	Premie (<28cm head circumference)	50/case
99071-02	Small (28-34cm head circumference)	50/case
99071-03	Medium (>34cm head circumference)	50/case



Customer Service: 1-800-345-6443 or 724-387-4000  
 Respironics Europe: +33-1-55-60-19-80  
 Respironics Asia Pacific: +81-3-5280-9611

WeeSpecs is a registered trademark of Respironics, Inc. People. Products. Programs. is a servicemark and Power Programs is a trademark of Respironics, Inc. and its affiliates. ©2004 Respironics, Inc.

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

1023291 KB 10/4/04 G

